



Introducing ViziVZV™

The first rapid lateral Flow Assay for detection of varicella (chickenpox) and zoster (shingles) IgG antibody.





Fast

In just 20 minutes, get results.



Easy to Use¹


Fingerstick



Prepare the test



Get Results



Accurate

**Sensitivity
100%**

**Specificity
100%**

**Detects VZV IgG
antigens even at
low levels not
detectable by ELISA**



Efficient

**No staff training, expensive equipment,
packaging, shipping, post-visit follow-ups.**



Features

- Single use cassette
- Includes lancets, capillary tubes and chase buffer
- Store at room temperature
- 12 month shelf-life
- Sold in boxes of 30 cassettes

Trusted Science

- Proprietary truncated VZV glycoprotein E (gE) capture antigen
- Most abundant and immunogenic VZV envelope glycoprotein
 - Capable of eliciting both humoral and cell-mediated immunity

**To order call 888-626-2047 or
e-mail us at info@viroresearch.com**

**For clinical questions
VIRORESEARCH.COM**

¹ See Instructions for Use for detailed directions.

