



Introducing ViziHSV™

The first rapid HSV1&2 Lateral Flow Assay for detection and differentiation of Herpes Simplex Virus types 1 and 2 (HSV1&2) IgG antibody.





Fast

In just 20 minutes, get results.



Easy to Use¹


Fingerstick



Prepare the test



Get Results



Accurate

Sensitivity

99%

Specificity

100%

Detects HSV IgG
antigens even at
low levels not
detectable by ELISA



Efficient

No staff training, expensive equipment,
packaging, shipping, post-visit follow-ups.



Features

- Single use cassette
- Includes lancets, capillary tubes and chase buffer
- Store at room temperature
- 12 month shelf-life
- Sold in boxes of 30 cassettes

Trusted Science

- Proprietary truncated HSV glycoprotein G (gG) capture antigen
- Most abundant and immunogenic VZV envelope glycoprotein
 - Capable of eliciting both humoral and cell-mediated immunity

¹ See Instructions for Use for detailed directions.

To order call 888-626-2047 or
e-mail us at info@viroresearch.com
For clinical questions
VIRORESEARCH.COM

For Investigational or Research Use Only



© 2023 Viro Research, LLC. All rights reserved in all countries.
Research and ViziVZV are trademarks of Viro Research, LLC.

04/2024